

# Owner's Manual

## 1000F Freestanding Trapeze



### ***Bariatric Freestanding Trapeze by ConvaQuip***

***Specializing in Bariatric Patient Aids***

P.O. Box 3417 Abilene, TX 79604 U.S.A  
Toll 800/637-8436 \* Fax 325-677-7217

# ConvaQuip Ind., Inc.

## PLEASE READ CAREFULLY!

- Do not install this trapeze while patient is in the bed !!
- Once trapeze has been assembled, be sure leg levelers are set to level the trapeze so it does not rock back and forth.
- Slide trapeze legs under the head end of the bed making sure the trapeze legs do not interfere with the bed frame and any electrical or manual/mechanical moving parts.
- Position trapeze frame in center of the bed.
- Position boom in center locator position so it is centered over the bed.
- Adjust boom length using slider on top of boom and secure with snap button.
- Adjust chain length to the desired position, making sure quick link is secured and chain can not disengage from handle.
- Boom can swivel 0 to 180 degrees. Be sure whatever position you choose, that the position pin is seated into slot so boom will not swivel while in use by the patient.

## Caution!

### Note:

**Position 0 and 180 Degrees are primarily used to get boom out of the way when entering and exiting the bed.**

**The 45 degree settings are commonly used for when the patient needs assistance when rolling over and not for entering and exiting the bed.**

**The 90 degree setting is commonly used to help the patient lift up or assist in sitting upright.**



## ConvaQuip Ind., Inc.

P.O. Box 3417  
Abilene, TX 79604  
Toll 800-637-8436  
Fax 325-677-7217

*Thank you for choosing ConvaQuip !*

## MAINTENANCE AND CLEANING



Capacity 1000 lbs.

## Cleaning:

- \* For infection control, clean all parts with one part bleach and 9 parts water. Rinse and allow to air dry or wipe with a clean, dry cloth.
- \* Pail may be steam cleaned or use Beaucoup disinfectant (made by EcoLab)



## ConvaQuip Ind., Inc.

P.O. Box 3417  
Abilene, TX 79604  
Toll 800-637-8436  
Fax 325-677-7217



# ConvaQuip Ind., Inc.

## Assembly Instructions for Model 1000F Freestanding Trapeze

Toll 800-637-8436 Fax 325-677-7217



Tools needed: Rubber Mallet  
3/4" Wrench & Drill



Parts: 2 Wheels, 2 Bolt Axles  
4 Levelers, 4 Self Tapping Screws  
and Socket Cap Screw



Lay base down as pictured  
and bolt on both wheels  
securely.



Screw each leveler all the  
way in. Do not over tighten.



Set base in upright position.



Align button on each leg and  
insert into top opening.



Make sure button secures  
in place.



Insert M-Shaped part into  
the 3 slots. Adjust all  
levelers so they touch floor.



Use rubber mallet to tap  
down until it bottoms out.



Using drill and socket cap  
screw provided, screw in the  
4 self tapping screws.



Insert Boom into slot on top  
of the M-Shaped piece.  
Make sure it bottoms out.



You can change positions of  
the boom by lifting pin



Unscrew the Quick-Link and  
hook over the Boom slider,  
then screw Quick-link closed



Trapeze is designed to tilt and  
roll. It is recommended to re-  
move the boom prior to moving.



Trapeze shown in most  
compact configuration



Easiest configuration for  
transporting from room  
to room.

# ConvaQuip Ind., Inc.

## LIMITED WARRANTY

This limited warranty is extended to the original purchaser/user and cannot be transferred. The trapeze frame is warranted for Lifetime against defects in material and workmanship. All other components are warranted for 1 year against defects in material and workmanship

Defective parts will be repaired or replaced by the company or its appointed agent. The warranty does not cover damage caused by misuse or negligence, nor does it cover defects or damages by use of unauthorized parts or service by an unauthorized person.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state.

Do not return merchandise without prior authorization. All correspondence relating to this warranty must indicate Model number and Serial number and be addressed to:

**ConvaQuip Ind., Inc.**  
**P.O. Box 3417**  
**Abilene, TX 79604**

## NOTICE

**FREIGHT DAMAGE** -The carrier who delivers merchandise to your door is responsible for loss and damages. Acceptance of the shipment by you is acknowledgment that all articles delivered were in good condition and properly packed. Therefore, all claims for loss or damage must be filed immediately with the freight carrier. Then notify us. We will mark our records accordingly. After settling claim, damaged merchandise will be picked up by the freight carrier and returned. Should you need assistance with the claim, call our customer service at 800-637-8436

**RETURN GOODS** - Within 60 days from date of purchase returns will be subject to a handling and re-stocking charge plus any additional costs resulting from damage or use. After 60 days from date of purchase: ConvaQuip Inc. reserves the right to not accept return goods for any reason. If accepted, the amount of credit will be at the discretion of ConvaQuip Inc. Returned merchandise without prior authorization will not be accepted.



**ConvaQuip Ind., Inc.**

P.O. Box 3417  
Abilene, TX 79604  
Toll 800-637-8436  
Fax 325-677-7217