



User Guide 8900

8900 USER INSTRUCTIONS



1. Place the commode on its side and find the 2 mounting brackets on the bottom of the seat.
2. Guide the tubes of the backrest through two of the mounting brackets of the seat.
3. Depress the push-button on the back tubes to allow them to slide into the seat bracket and push the tubes in until both push-buttons engage into the bracket holes.
4. Adjust the height of the legs by depressing the push-button and sliding the leg up or down to the desired height. Make sure both legs are fully engaged and adjusted to equal heights.
5. **RE-CHECK ALL PUSH BUTTONS TO MAKE SURE THEY ARE FULLY ENGAGED BEFORE USING.**

How to Use

To Drop Arm – Push gray button located on the side of the arm (see above)
 Pull arm up until click is heard and arm is secure before bearing weight

Over a toilet - Remove bucket and lid from the Bedside Commode.
 Lift up the toilet seat and lid.
 Place Bedside Commode over the toilet.

As the toilet - Make sure the bucket is placed under the commode lid.
 Remove the bucket lid.
 Make sure toilet seat is lowered to a comfortable height

Care and Maintenance

Check regularly to make sure the back and legs are securely set.
 Clean regularly with a toilet sanitizer or cleanser.
 Discard bucket contents immediately after use.
 Clean bucket immediately after use.

*******PATIENT WEIGHT NOT TO EXCEED 300LBS.*******

LIFETIME LIMITED WARRANTY

Nova Ortho-Med, Inc. warrants this product to be free from defects in material or workmanship for a period of lifetime from date of purchase. If within the specified warranty period, this product is proven to be defective, it will be repaired or replaced at NOVA O.M.'s option. This warranty does not cover shipping costs incurred in replacement or repair of this product. NOVA O.M. stands firm on our commitment to providing the highest quality products, as well as exceptional service. This product was built to precise standards and thoroughly inspected prior to shipment. This Lifetime Limited Warranty represents our confidence in the materials and workmanship of our products. Do not return any products to our office without prior consent. This warranty does not cover problems due to user negligence or misuse of the product. This warranty shall not extend to non-durable components, including but not limited to rubber tips, handgrips, casters, and other parts subject to normal "wear and tear" THIS FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS WARRANTIES, IMPLIED WARRANTIES, IF ANY, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NOVA O.M. SHALL NOT BE LIABLE FOR ANY SEQUENTIAL OR INCIDENTAL DAMAGES WHATSOEVER.