
Mini Bike EXERCISER

Components List

- ▶ 1 - Main Body
- ▶ 2 - Front Stabilizer
- ▶ 3 - Rear Stabilizer
- ▶ 4 - Right Pedal
- ▶ 5 - Left Pedal
- ▶ 6 - Monitor
- ▶ 7 - Tension Knob
- ▶ 8 - Crank Arm
- ▶ 9 - Non-skid Mat
- ▶ 10 - Wrench
- ▶ 11 - Hex Bolts (4) — Not Shown
- ▶ 12 - Washers (4) — Not Shown

PLEASE READ PRIOR TO USE

Safety Summary

- ▶ Consult your physical therapist or health care professional before starting any exercise program.
- ▶ Establish your own exercise routine, (frequency, intensity and time) that best suits your particular needs and condition.
- ▶ If you experience pain or tightness in your chest, an irregular heartbeat, shortness of breath, feel faint or have any discomfort while exercising, STOP immediately! Consult your physician or health care professional before continuing any exercise routine.

WARNING: Do not stand on the Mini Bike.

CAUTION: The metal components on the tension knob may become hot when in use.

Assembly Instructions

NOTE: Tighten all hex bolts (11) and pedals (4,5) with the included wrench (10).

- ▶ 1. Attach the front stabilizer (2) to the main body (1) with two hex bolts (11) and two washers (12), Fig. 1.
- ▶ 2. Attach the rear stabilizer (3) to the main body (1) with two hex bolts (11) and two washers (12), Fig. 1.
- ▶ 3. The right pedal (4), marked "R" on the end of the pedal, connects to the crank arm (8), Fig. 2. **NOTE: The right pedal should be threaded on clockwise. Do not over tighten as it could strip the threads.**
- ▶ 4. The left pedal (5), marked "L" on the end of the pedal, connects to the crank arm (8), Fig. 3. **NOTE: The left pedal should be threaded on counter-clockwise. Do not over tighten as it could strip the threads.**
- ▶ 5. Make sure all bolts are tightened securely before use.

WARNING: Periodically check assembly for tightness. Do not use if any parts appear to be loose or damaged as personal injury may result.

NOTE: If any parts are missing do not return to place of purchase. Contact Customer Care at 877-505-6190.

Fig. 1

Fig. 2

Fig. 3

Battery Installation / Replacement

- ▶ 1. Remove the monitor (6) from the top of the main body (1), Fig. 4.
- ▶ 2. Insert or replace one AAA battery in the compartment, matching the indicated polarity symbols.
- ▶ 3. Replace the monitor back in the top of the main body.

NOTE: When removing monitor from the top of the main body the wire connecting the monitor to the Mini Bike may become disconnected. Before replacing the monitor reconnect the wire to the Mini Bike, Fig. 5.

Removing the battery will automatically reset all the stored readings.

Fig. 4

Fig. 5

Exercise Monitor Modes

The red Mode button on the monitor (6) allows the user to switch between the various modes of the exercise monitor.

- ▶ Scan - Automatically scans through each mode in sequence.
- ▶ Time - Counts upwards to 99:59 from zero in one second increments.
- ▶ Calories - Accumulates an estimated total calorie consumption during exercise.
- ▶ Count - Accumulates the total number of revolutions during different sessions of your workout.
- ▶ Total Count - Accumulates the total number of revolutions during different sessions of your workout. **Count cannot be reset unless you replace the batteries.**

NOTE: To reset all readings to "0", press and hold the red Mode button for 4 seconds. This will reset all modes except for the Total Count. To reset the Total Count, remove the battery.

Monitor operates with one (1) AAA battery included. Monitor will automatically shut-off after 4 minutes.

Instructions For Use

The Mini Bike is ideal for individuals who are unable to sit or ride on a regular bicycle, or for those limited to wheelchair use. The Mini Bike can be used on the floor for leg exercises or on a table top for arm exercises.


Floor Use:

- ▶ 1. Place the non-skid mat (9) between the Mini Bike and the floor to help prevent slipping.
- ▶ 2. Adjust the tension knob (7) to the desired level of resistance. Turn the tension control knob clockwise to increase resistance and counter-clockwise to reduce resistance.
- ▶ 3. Once seated, adjust the Mini Bike to a comfortable position in front of you.
- ▶ 4. Place each foot on the pedals in the straps and begin pedaling.
- ▶ 5. Pedal forward or backward to work different muscles groups.

Table Top Use:

- ▶ 1. Place the non-skid mat (9) between the Mini Bike and the table to help prevent slipping.
- ▶ 2. Adjust the tension knob (7) to the desired level of resistance. Turn the tension control knob clockwise to increase resistance and counter-clockwise to reduce resistance.
- ▶ 3. Sit with the Mini Bike comfortably positioned in front of you. Keep your back straight.
- ▶ 4. Grasp each pedal with your hands and begin pedaling, trying not to lean forward.
- ▶ 5. Pedal forward or backward to work different muscles groups.

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©2009 Duro-Med Industries