

# ALPS™ Skin Lotion

# SL108

## Prosthetic Skin Lotion for Sensitive Skin



ALPS™ 100% Silicone Skin Lotion helps to soothe dry, sensitive skin. It is hypoallergenic, with no perfumes or dyes, and does not leave an oily feeling after use. Can be used sparingly to relieve dry, chapped skin conditions.

ALPS™ Skin Lotion is available in 4 oz. (SL108) and 32 oz. (SL108-32) size bottles.

Recommended for use with:

- ✓ Gel (thermoplastic) - All ALPS™ Gel Liners
- ✓ Silicone - CLEARPRO
- ✓ Urethane

# ALPS™ Prosthetic Lubricant

# PA604

## 70% Isopropyl Alcohol Prosthetic Lubricant for Liners



ALPS™ Prosthetic Lubricant is used to help in donning silicone prosthetic liners. It allows the patient to don the liner or sheath easily while eliminating the sticky feel.

ALPS™ Prosthetic Lubricant is available in a 16 oz. bottle.

Recommended for use with:

- ✓ Gel (thermoplastic) - All ALPS™ Gel Liners
- ✓ Silicone - CLEARPRO
- ✓ Urethane

# ALPS™ Prosthetic Cleanser

# PD595

## Prosthetic Cleanser for Sensitive Skin



ALPS™ Prosthetic Cleanser is a non-alkaline, soap-free cleanser which can be used for either hand or machine washing of prosthetic liners or sleeves. Patients with sensitive skin may also find that this soothing cleansing agent is ideal for washing skin areas that come into contact with the prosthesis.

ALPS™ Prosthetic Cleanser is available in a 16 oz. bottle.

Recommended for use with:

- ✓ Gel (thermoplastic) - All ALPS™ Gel Liners
- ✓ Silicone - CLEARPRO
- ✓ Urethane



Tel: 727.528.8566 Tel: 800.574.5426 [www.easyliner.com](http://www.easyliner.com)  
[info@easyliner.com](mailto:info@easyliner.com)