

Walkabout Four-Wheel Imperial Rollator

The Lumex Walkabout Four-Wheel Imperial (Bariatric) is ideal for larger individuals. Similar to the Walkabout Lite, this rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 20 lbs and supports a maximum weight capacity of 500 lb - **EVENLY DISTRIBUTED**. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Other features include:

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded 20" wide seat for increased seating surface and added comfort
- Removable padded backbar
- Adjustable handle height to accommodate broad range of user heights
- Durable, heavy duty 8" wheel for indoor and outdoor use
- Basket included that can be placed underneath seat for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0149 & E1056
- Includes easy to read, retail-friendly user manual
- Lifetime Limited Warranty

And now with added rear reflectors that increase visibility at night!

Straight Backbar

| | | |
|----------------|----------|------|
| RJ4400K | Black | 1 ea |
| RJ4400R | Burgundy | 1 ea |

Imperial Hemi

| | | |
|----------------|----------|------|
| RJ4402R | Burgundy | 1 ea |
|----------------|----------|------|

Contoured Backbar

| | | |
|----------------|------------|------|
| RJ4405B | Royal Blue | 1 ea |
| RJ4405K | Black | 1 ea |
| RJ4405R | Burgundy | 1 ea |

| Specifications | RJ4400 | RJ4405 |
|-----------------------------------|--------|--------|
| Width Inside Handgrips | 24" | 24" |
| Overall Width | 29 ½" | 29 ½" |
| Overall Depth | 27 ½" | 30" |
| Seat Width | 20" | 20" |
| Seat Depth | 13 ½" | 13 ½" |
| Seat Height | 23 ½" | 23 ½" |
| Caster Diameter | 8" | 8" |
| Minimum Handle Height | 34" | 34" |
| Maximum Handle Height | 39" | 39" |
| Weight Capacity | | |
| - EVENLY DISTRIBUTED | 500 lb | 500 lb |
| Product Weight | 20 lbs | 20 lbs |



Burgundy model shown



Color options

Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.