

Pelvic Stabilizer



HCPCS
E1399

Item #	For Use With	UOM
CE 1050	All Wenzelite Safety Rollers	1/bx

- The pelvic stabilizer is adjustable in width and features a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and, therefore, cannot stay centered in the walker.

Vertical Handgrips



HCPCS
E1399

Item #	For Use With	UOM
CE 1045	All Wenzelite Posterior Safety Rollers	1/bx

- Patients who cannot grasp a standard handlebar use vertical handgrips. The vertical handgrips extend, pronate, supinate and have a 45-degree angle. Each side is independently adjusted for the most comfortable hand placement. When both vertical handgrips are placed in the prone position, they can be used as a horizontal handlebar.

Positioning Bar



HCPCS
E1399

Item #	For Use With	UOM
CE 1052	PE TYKE 1200	1/bx
CE 1053	PE 1200	1/bx
CE 1054	CE 1200 BK	1/bx

- The positioning bar is mounted on the handlebar of Safety Rollers and Gait Trainers and is easily moved for adjustable handlebar depth. It can be used with or without the pelvic stabilizer.

Seat Harness



HCPCS
E1399

Item #	Product Size	UOM
CE 1070S	Small	1/bx
CE 1070L	Large	1/bx

- The Seat Harness provides cushioned comfort. It moves with user's gait pattern allowing the user to shift weight from one leg to the other. The straps can be lengthened, shortened and positioned towards the front or back on the walker frame. It assists in weight-bearing, properly positions user's pelvis and can be used as a safety sling.

For additional specifications, Page 60.