


# POSITIONING AND BED CARE



## BED WEDGE




### FEATURES

- May be used in a variety of positioning applications.
- Solace cover.

| BED WEDGE    |                 |
|--------------|-----------------|
| MODEL NUMBER | L X W X H       |
| <b>88A</b>   | 22" x 20" x 18" |
| <b>88B</b>   | 25" x 20" x 14" |
| <b>88C</b>   | 26" x 20" x 10" |
| <b>88D</b>   | 28" x 20" x 6"  |

## BODY BOLSTER



### FEATURES

- Full length positioning device used in therapy or in bed.
- Solace cover.

| BODY BOLSTER |                 |
|--------------|-----------------|
| MODEL NUMBER | L X W X H       |
| <b>87</b>    | 32" x 12" x 12" |

## SOFT RAIL



### FEATURES

- Firm foam provides safety from bed rail.
- Easily cleaned solace fabric cover.
- Can be mounted anywhere on bed.

| SOFT RAIL    |               |
|--------------|---------------|
| MODEL NUMBER | L X W X H     |
| <b>92</b>    | 32" x 7" x 8" |

SOLD AS SINGLE UNIT

## SIT 'N PLACE



STANDARD



SUPER

### FEATURES

- Connects to chair with nylon strapping system and side release buckle.
- Solace cover.
- Larger sized Super Sit 'n Place is available.

| SIT'N PLACE STANDARD |            | SIT'N PLACE SUPER |            |
|----------------------|------------|-------------------|------------|
| MODEL NUMBER         | DIMENSIONS | MODEL NUMBER      | DIMENSIONS |
| <b>100SP</b>         | STANDARD   | <b>200SSP</b>     | SUPER      |

## SAFE & SECURE



### FEATURES

- Steel clamps and screws secure Safe & Secure to bed rail with no straps.
- Solace cover.
- Full rail or half rail lengths are available.

| SAFE & SECURE FULL RAIL |                | SAFE & SECURE HALF RAIL |                |
|-------------------------|----------------|-------------------------|----------------|
| MODEL NUMBER            | L X W X H      | MODEL NUMBER            | L X W X H      |
| <b>83A</b>              | 64" x 2" x 14" | <b>83B</b>              | 32" x 2" x 14" |

# POSITIONING AND BED CARE



## ABDUCTION PILLOW



### FEATURES

- Enables users to maintain proper leg positioning while recovering from hip replacements or fractures.
- Double padded straps for a secure fit.
- Solace cover.

### ABDUCTION PILLOW

| MODEL NUMBER | L X W X H      |
|--------------|----------------|
| <b>73A</b>   | 14" x 14" x 7" |
| <b>73B</b>   | 18" x 18" x 7" |

## LEG ABDUCTOR WEDGE



### FEATURES

- Enables users to maintain proper leg positioning while recovering from hip replacements or fractures.
- Connects easily to user with soft hook and loop straps.
- Solace cover.

### LEG ABDUCTOR WEDGE

| MODEL NUMBER | L X W X H    |
|--------------|--------------|
| <b>72</b>    | 7" x 9" x 6" |

## HIP BOLSTER



### FEATURES

- Protects user from side of chair.
- Solace cover.
- Strapping system secures around wheelchair side.

### HIP BOLSTER

| MODEL NUMBER | L X W X H     |
|--------------|---------------|
| <b>81A</b>   | 14" x 2" x 6" |

### HIP BOLSTER

| MODEL NUMBER | L X W X H     |
|--------------|---------------|
| <b>81B</b>   | 18" x 3" x 6" |

## GAIT BELT

HCPCS CODING: E0978



FAST ATTACH:  
METAL CONNECTION



SIDE RELEASE:  
PLASTIC CONNECTION

### FEATURES

- Heavy-duty 2" cotton webbing.
- Available with a plastic side release connection or fast attach metal connection

### GAIT BELT FAST ATTACH METAL

| MODEL NUMBER | L X W    |
|--------------|----------|
| <b>97A</b>   | 49" x 2" |
| <b>97B</b>   | 60" x 2" |
| <b>97C</b>   | 73" x 2" |

### GAIT BELT SIDE RELEASE PLASTIC

| MODEL NUMBER | L X W    |
|--------------|----------|
| <b>96A</b>   | 49" x 2" |
| <b>96B</b>   | 60" x 2" |
| <b>96C</b>   | 73" x 2" |