

Everest & Jennings®  
**DURA-gel™**  
**SPP**

The Dura-Gel™ SPP (Skin Protection & Positioning) wheelchair cushion reduces surface pressure through the use of non-aqueous Dura-Gel™ placed in three pressure zones. The cushion uses a high-density base foam in combination with a visco-elastic topper for added comfort and also positions the user with a moderate slope from front to back. Cushion also features medial thigh positioning as well as lateral positioning to provided greater stability for the user.



- ❑ Provides long-term skin protection and positioning for higher risk patients
- ❑ Reduces tissue interface pressures at critical points
- ❑ Features a specialized Dura-Gel™ center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ❑ Side wedges provide lateral stability to the lower extremities, trochanter and pelvis. Center pommel provides medial thigh support
- ❑ The soft, stretchable lightweight cover reduces shear/friction and is fluid resistant and has a strap at the back to secure the cushion firmly to the chair
- ❑ Accommodates a wide range of body shapes, sizes and weights
- ❑ 350 lb maximum weight capacity
- ❑ Lightweight
- ❑ 2 year limited warranty



**Latex Free**

*Designed to meet the requirements of HCPCS Code: E2607*

- 8920166** 16" x 16" x 4 ½"
- 8920186** 18" x 16" x 4 ½"
- 8920188** 18" x 18" x 4 ½"
- 8920206** 20" x 16" x 4 ½"
- 8920208** 20" x 18" x 4 ½"

*Designed to meet the requirements of HCPCS Code: E2608*

- 8920228** 22" x 18" x 4 ½"
- 8920248** 24" x 18" x 4 ½"

