

Forearm Supports for KAYE Posture Walkers

Models FSS, FSM, FSH

Description and Use

Kaye Forearm Supports are available for Kaye Posture Control, PostureRest and Red walkers. The handgrips are adjustable from horizontal to vertical to accommodate supination or pronation of the forearm. The forearm pads adjust in height as well as forward-backward to accommodate shoulder and arm positioning. These pads can be angled to change the alignment of the shoulder to the elbow.

Features and Specifications

Forearm Support Model Number	Height from top of walker handle to top of forearm cushion.	Fits Walker Models – all wheel configurations
FSS	4" to 7 ³ / ₄ "	W1/2B, W1/2 BH, W1B, W1BH, W2B, R2B, W3B, R3B
FSM	6 ¹ / ₂ " to 10"	W4B, R4B, W5C, W5CH
FS2H	4" to 7 ¹ / ₂ "	W2H
FS3H	6" to 9"	W3H
FS4H	5" to 9"	W4H



FSS shown on W2BS walker



FSA on RF-B walker



Model FS2H on KAYE PostureRest Walker model W2H.



Model FSS

W1/2BS walker