

# THE SUPPORT PRO 7 SERIES

GENERAL USE CUSHION

HCPCS CODING

E2602




WOUND STAGE 3

## PRODUCT FEATURES

Kwik Strap® loops around vertical canes to securely attach cushion to wheelchair during transfers.

Non-contoured seating surface allows for easy transfers.


Coccyx relief area unloads bony prominences.



Waterfall front eliminates pressure points in leg trough.

750 lb. weight capacity foam provides increased comfort and reduces potential high-pressure areas.

Optional QuadraGel® pack insert.



SEE CHART FOR DIMENSIONS

### COVER STYLE



**SOLACE FABRIC:** is designed for pressure relief and comfort. Solace is fluid-proof, anti-microbial and soft to the touch. The multi-dimensional stretch design alleviates pressure on bony prominences and provides optimal skin protection.



### QUADRAGEL® PACK

The QuadraGel® pack regulates body temperature and redistributes pressure for comfort.



#### SUPPORT PRO 7 SERIES WITHOUT QUADRAGEL

SOLACE	W X L	FRONT	REAR	WT. LIMIT
<b>53-2218-B</b>	22 x 18	3.5"	3.5"	750 lbs
<b>53-2220-B</b>	22 x 20	3.5"	3.5"	750 lbs
<b>53-2418-B</b>	24 x 18	3.5"	3.5"	750 lbs
<b>53-2420-B</b>	24 x 20	3.5"	3.5"	750 lbs
<b>53-2618-B</b>	26 x 18	3.5"	3.5"	750 lbs
<b>53-2620-B</b>	26 x 20	3.5"	3.5"	750 lbs
<b>53-2818-B</b>	28 x 18	3.5"	3.5"	750 lbs
<b>53-2820-B</b>	28 x 20	3.5"	3.5"	750 lbs
<b>53-3018-B</b>	30 x 18	3.5"	3.5"	750 lbs
<b>53-3020-B</b>	30 x 20	3.5"	3.5"	750 lbs

#### SUPPORT PRO 7 SERIES WITH QUADRAGEL

SOLACE	W X L	FRONT	REAR	WT. LIMIT
<b>53G-2218-B</b>	22 x 18	3.5"	3.5"	750 lbs
<b>53G-2220-B</b>	22 x 20	3.5"	3.5"	750 lbs
<b>53G-2418-B</b>	24 x 18	3.5"	3.5"	750 lbs
<b>53G-2420-B</b>	24 x 20	3.5"	3.5"	750 lbs
<b>53G-2618-B</b>	26 x 18	3.5"	3.5"	750 lbs
<b>53G-2620-B</b>	26 x 20	3.5"	3.5"	750 lbs
<b>53G-2818-B</b>	28 x 18	3.5"	3.5"	750 lbs
<b>53G-2820-B</b>	28 x 20	3.5"	3.5"	750 lbs
<b>53G-3018-B</b>	30 x 18	3.5"	3.5"	750 lbs
<b>53G-3020-B</b>	30 x 20	3.5"	3.5"	750 lbs